

#####  
#  
# Bruker BioSpin now supports  
# Windows XP Professional SP3 and  
# Windows Vista SP1 with TopSpin  
#  
#####

Microsoft has released the official service pack 1 for Windows Vista and service pack 3 for Windows XP.

After detailed testing and validation of both Windows service packs on different spectrometers with TopSpin software, no limitation of spectrometer handling and TopSpin use could be diagnosed.

Therefore, Bruker BioSpin now supports Windows Vista Business service pack 1 and Windows XP Professional service pack 3 with TopSpin for acquisition and for processing.

Please note that Windows Vista requires TopSpin 2.1.

For further information about software compatibility, please refer to the Software Support Matrix under the following Bruker Web address:  
[http://www.bruker-biospin.com/documentation\\_general.html](http://www.bruker-biospin.com/documentation_general.html)

-----  
Bruker BioSpin GmbH  
Silberstreifen 4  
76287 Rheinstetten  
Germany  
-----

To cancel your subscription to the Bruker NMR software newsletter, reply to this e-mail with the word UNSUBSCRIBE in the subject line.

To cancel your subscription to all Bruker BioSpin newsletters, reply to this e-mail with the word UNSUBSCRIBE\_ALL in the subject line.