

```
#####  
#  
#           Release of TopSpin 2.1. pl4           #  
#  
#####
```

The new TopSpin patch level 4 for version 2.1 (2.1pl4) is now available for:

Acquisition:

- Windows XP Professional
- Windows Vista (Business Edition)
- Red Hat Enterprise Linux WS 3
- Red Hat Enterprise Linux WS 4

Processing:

- Windows 2000
- Windows XP Professional
- Windows Vista (Business Edition)
- Red Hat Enterprise Linux WS 3
- Red Hat Enterprise Linux WS 4
- Red Hat Enterprise Linux WS 5

TopSpin 2.1pl4 comes as a full program DVD and contains the following program versions:

- TopSpin 2.1pl4
- TopSpin Plot Editor 4.5
- ICON-NMR 4.2.4
- NMR-SIM 5.0.2
- NMR-GUIDE 4.0
- AUTOLINK

The following programs are distributed as part of the TopSpin 2.1pl4 DVD, but they must be licensed separately:

- AMIX 3.8.7
- AUREMOL 2.3.1
- PERCH NMR TOOLS 2007-09-10 (only for Windows operating system)

TopSpin 2.1pl4 can be used with all AVANCE, AVANCE-II and AVANCE-III spectrometers. If you are not sure about upgrading possibilities or your hardware compatibility, please refer to the Release Letter of TopSpin 2.1pl4 for information about EC-levels or take a look at the Software Support Matrix for hardware components under the following Bruker web address:

[http://www.bruker-biospin.com/documentation\\_general.html](http://www.bruker-biospin.com/documentation_general.html)

TopSpin 2.1pl4 is a minor upgrade from TopSpin 2.1 with many improvements concerning bug fixes.

Please order the TopSpin 2.1pl4 DVD from your local Bruker office. Bruker BioSpin addresses can be found here:

[http://www.bruker-biospin.com/contact\\_us.html](http://www.bruker-biospin.com/contact_us.html)

It is recommended by Bruker BioSpin to install TopSpin 2.1pl4 in a parallel directory to your existing TopSpin 2.1 or a previous version. Please note that TopSpin 2.1pl4 contains the full program, so it is not necessary to install TopSpin 2.1 or a previous version first.

For full information of all changes, bug fixes and FAQs, please refer to the TopSpin 2.1pl4 Release Letter, which is available electronically on the TopSpin 2.1pl4 DVD (in pdf and rtf format) and on the Bruker website under the following address:

[http://www.bruker-biospin.com/documentation\\_general.html](http://www.bruker-biospin.com/documentation_general.html)

All bug fixes and FAQs can also be found in the Bruker Knowledge Base on the Bruker Website under the following address:

[http://www.bruker-biospin.com/bkb\\_access.html?&no\\_cache=1](http://www.bruker-biospin.com/bkb_access.html?&no_cache=1)

Please note that a user account and password is required to open this page. To login or register for a free account, please open the following web page:

<http://www.bruker-biospin.com/register.html>

If you are already connected to the internet and if you have access to the Bruker Knowledge Base, you can then search for all bug fixes in TopSpin 2.1pl4 or you can enter the item number in the appearing window.

Furthermore there are several important FAQs available in the Bruker Knowledge Base. In the Knowledge Base you have to select program "TopSpin" and status "FAQ" to enter the FAQ-area that will give you always the up-to-date list of all interesting TopSpin questions respectively answers.

Please note that for entering the NMR Software Knowledge Base you always need a website login.

Please note that the most recent version of the Software Release Letter for TopSpin 2.1 pl4 can be found on the Bruker website at the following link:

[http://www.bruker-biospin.com/documentation\\_general.html](http://www.bruker-biospin.com/documentation_general.html)